

## Product Description



One-time Used 134.2Khz FDX-B 2.12x12mm EM4305 RFID Pet Tracker Glass Tag Animal Tracking Pet Microchip For Dogs Fish

RFID Glass Tag is a flexible and reliable solution for pet, livestock and fish identification, pigeon race and identification of laboratory animals. It is a flexible solution for secure animal identification that withstands most chemicals and harsh environment circumstances .

Features:

- ◆ Especially designed for animal identification
- ◆ Bio-compatible glass encapsulation
- ◆ Small size with long read range required
- ◆ Available anti-collision in LF and HF
- ◆ Unlimited resistance to water and chemical absorption
- ◆ High stability over fluctuating temperatures

## Specification



S  
p  
e  
c  
i  
f  
i  
c  
a  
t  
i  
o  
n

R  
F  
b  
D  
A  
n  
t  
m  
a  
m  
G  
l  
a  
s  
s  
T  
u  
b  
e  
T  
a  
g  
M  
i  
o  
g  
r  
i  
a  
s  
E  
M  
#  
D  
p  
N  
T  
A  
G  
2  
1  
6

2  
i  
1  
2  
\*  
1  
2  
m  
m  
/  
1  
.  
2  
5  
\*  
7  
m  
m  
/  
1  
.  
4  
\*  
8  
m  
m

A  
p  
p  
h  
à  
c  
h  
d  
è  
n  
h  
i  
f  
i  
c  
a  
t  
i  
o  
n  
o  
r  
t  
r  
a  
c  
k  
i  
n  
g  
W  
b  
fC  
k  
b  
fC  
g  
t  
e  
m  
p  
e  
r  
a  
t  
u  
r  
e

Part of our HF chips, for more RFID chips (LF, HF, UHF, and dual frequency, please contact us freely)

Protocol ISO/IEC 14443A:

1. MIFARE Classic® 1K, MIFARE Classic® EV1 1K, MIFARE Classic® 4K

MIFARE and MIFARE Classic are registered trademarks of NXP B.V. and are used under license.

2. MIFARE Plus® S 1K, MIFARE Plus® S 1K SE, MIFARE Plus® S 2K / S 4K, MIFARE Plus® X 2K / X 4K,  
MIFARE

Plus® EV1 2K / 4K

MIFARE and MIFARE Plus are registered trademarks of NXP B.V. and are used under license.

3. MIFARE® DESFire®

MIFARE® DESFire® 2K / EV1 2K / EV2 2K

MIFARE® DESFire® 4K / EV1 4K / EV2 4K

MIFARE® DESFire® 8K / EV1 8K / EV2 8K

MIFARE DESFire are registered trademarks of NXP B.V. and are used under license.

4. NFC Forum Type 2:

1) NTAG® 203 (144 bytes), NTAG® 213 (144bytes), NTAG® 215 (504 bytes), NTAG® 216(888 bytes),  
NTAG® 210 (48 bytes), NTAG® 212 (128

bytes) NTAG® are registered trademarks of NXP B.V. and are used under license.

2) MIFARE Ultralight® (48 bytes) MIFARE Ultralight® EV1 (48 bytes/128 bytes) MIFARE Ultralight® C(148  
bytes)

MIFARE and MIFARE Ultralight are registered trademarks of NXP B.V. and are used under license.

Protocol ISO 15693/ISO 18000-3:

ICODE® SLIX, ICODE® SLIX-S, ICODE® SLIX-L, ICODE® SLIX 2

ICODE® are registered trademarks of NXP B.V. and are used under license.

Details Images

# How Does Animal Microchip Work?



1. Scan the pet to verify whether he or she has a microchip that has been previously implanted.



2. Scan the microchip pack to verify that the microchip is functional, and that the ID code matches the accompanying barcode labels.



3. Implant the microchip, then rescan the pet to verify that the microchip can be localized.



4. Do not forget to record the tracking number of the microchip.

## Implantation Sites for Microchips



Microchips are most often implanted in animals as outlined in the diagrams above. When scanning dogs and cats, begin in the neck area where the microchip is most likely to be located, but be sure to scan slowly and patiently. Scan repeatedly and over the entire body to ensure whether or not a microchip can be localized.

# WHY YOU CHOOSE US AS YOUR PARTNER



## NO. 1 MATERIAL

Material is medical grade bio-glass with working life as long as 20 years

## NO. 2 STERILIZATION

Sterilized by EO for 48 hours to protect animal from being infected.

Safe enough to be injected into animal body.





# NO. 3 CHIP

Imported original chip. Long reading range and long working life




ISO11784/5, FDX-B, Chip Type: EM4305

ISO14443A/15693, Chip Type: NFC 216, etc

# NO. 4 VARIOUS SIZES






**NO. 5**  **VARIOUS MODELS AND COLORS**  
SYRINGES FOR ANIMAL ID TAGS



Application

# APPLICATION



Fish tracking



Big Animal ID



Pet ID



Small animal ID



Bird ID



Packing & Delivery



# PACKING



## Our Service

- 1, Any inquires will be replied within 24 hours
- 2, Professional manufacturer and supplier, Welcome to visit our website and our factory
- 3, OEM/ODM Available
- 4, High quality, fashin desing, reasonable & competitive price, fast lead time
- 5, After-sale Service :
  - 1), All products will have been strictly quality checked in house before packing
  - 2), All products will be well packed before shipping
  - 3), All our products have 2-3 year warranty if the damage is not caused by human
- 6, Faster delivery : Around 1~ 5 days for sample order, 7~30 days for bulk order
- 7, Payment : You can pay for the order via : T/T, Western Union, Paypal, Ali trade assurance
- 8, Shipping: We have strong cooperation with DHL, FEDEX, TNT, UPS, EMS, Forwarder by SEA and By AIR, You also can choose your own shipping forwarder.

### Quality Warranty

- Warranty Service will be honored if the damage is not caused by human, ACM Goldbridge provides 2 years warranty for relative products.
- On the contrary, ACM Goldbridge will cahрге extra if repair.
- More information, please browse our service center.

## Our Company

ACM Shenzhen Goldbridge industrial company

As a one of leading manufacturer and exporter of RFID and access control and electric door locks products in China , We have been in this field for 20 years already . We have rich production and exportation experience on RFID access control locks products . Our strength products are : RFID card , RFID keyfob , RFID wristband , RFID tag and various RFID reader, UHF long range reader, electric door locks, . We are also the access control solution provider .



## Contact us

Please send us inquiry and order now to have your FREE SAMPLE in 3 days!!!

We are professional and make sure can offer best solution of comprehensive security rfid card access control products supplier with 20 years experience,

So, Contact us and send your inquiry details in the below now!

We are looking forward to forming successful business relationships with new clients around the world in the future !



**SHENZHEN GOLDBRIDGE  
INDUSTRIAL CO., LTD**

Address:

Building 405, Shangbu Industrial  
Area, Futian District, Shenzhen China,  
518031

Tel: +0086-755-83767678 Fax:  
+0086-755-88391348

Contact Person: Lily Jiang

Mob/whatsapp/wechat:  
+0086-13554918707

Skype: lily-jiang1206

E-mail: sales(at)goldbridgesz.com

Website:

<http://www.goldbridgesz.com>